



ONBOARDING JOURNEY - SIMPLIFIED

Register with Brandsync

- 1. Register with
 Brandsync and create
 your company's
 account
- 2. Obtain a GLN for every location / warehouse / store that deals with medical products

To get a GLN you need to Contact: tatmeen@gs1-uae.org

Upon registration
Tatmeen will receive
your company details
and your products'
details and will send you
an invitation to register
for the testing system.

Register with Tatmeen Staging Environment

Upon receiving the invitation email, you need to:

- 1. Register your username and password on Tatmeen Staging System. this will be the single point of contact for your organization (SPOC)
- 2. Review your company details and product details

For Mobile & Portal Users

Once you register and create SPOC account you need to:

- 1. Create a "Mobile /
 Portal" user. This
 user will be able to
 access Tatmeen
 portal (website) and
 Mobile application
- 2. Familiarize yourself with the system by reading the manuals published on Tatmeen.ae

If you are a MAH or
Licensed Agent then you
will see your products on
the portal, you need to
ensure that the product
data is correct

Step 4 is applicable only for partners who plan to connect their system with Tatmeen using B2B connectivity

For Business
To Business
(B2B)

For B2B connection, you need to:

- 1. Create a "Service
 User" account to be
 able to use the
 Tatmeen API
 (Application
 Programming
 Interface)
- 2. Connect your stock management system to Tatmeen using API
- 3. Perform Test case scenarios as explained in the technical guidelines to confirm successful connectivity

Connect to Live System

Towards the Go Live date. Tatmeen will send the SPOC a form to check the method of connection.
All partners need to repeat steps 2 and 3 on the live system.

Partners who decided to connect via B2B will need to provide proof of successful test cases as published on Tatmeen.ae before they get the live system access. And then they need to repeat steps 2-4 and update their system accordingly.



GO LIVE

13th DEC 22

